

Nostril Retainer

This product prevents the deformation of the nasal foramen after the surgery for cleft lip or palate, rhinoplasty, or augmentation rhinoplasty without imposing a burden on the patient



UNOSTRIL

Nostril Retainer



EASY USE WITH COMFORT



Nostril Retainer UNOSTRIL

- This product prevents the deformation of the nasal foramen after the surgery for cleft lip or palate, rhinoplasty, or augmentation rhinoplasty without imposing a burden on the patient.
- It is available in 8 sizes each so that it can be used in a size appropriate to each patient.
- It is made of soft, less irritative silicone.

Non-Sterile

The Nasal retainer -also known as a nostril retainer or a nasal stent- can be used separately or as part of the NAM-treatment (Nasoalveolar Molding) for infants born with cleft lip syndrome.

Prior to surgery, the device helps to correct deformities of the nasal septum and nostrils to a more centered shape. After surgery, UNOSTRIL is used to sustain the results of the surgical correction and further shape the nasal dome

Integral airway allows post-operative nasal breathing

The UNOSTRIL Nasal Retainer shapes the nasal dome in a natural centered position and creates good airflow at the same time. The curved design enhances the symmetrical shape of the nostrils. The device has two tubes that are interconnected with a 'connection bridge' for safety.

We use a special biocompatible silicone that combines optimal user comfort (soft touch) with enough resilience to shape the nose.

Contact Us

Genco Tibbi Cihazlar San. Tic. Ltd. Sti

1571 Sok. No:22 Kat:2
Cinarli/Izmir/TURKEY
T: +90 232 4621649
F: +90 232 4353116
info@gencotibbi.com

www.gencotibbi.com

TABLE OF SIZES:

PRODUCT NO	INTRODUCTION		BARCODE NO
NST01	NOSTRIL RETAINER NO: 1		8680614811081
NST02	NOSTRIL RETAINER NO: 2		8680614811098
NST03	NOSTRIL RETAINER NO: 3		8680614811104
NST04	NOSTRIL RETAINER NO: 4		8680614811111
NST05	NOSTRIL RETAINER NO: 5		8680614811128
NST06	NOSTRIL RETAINER NO: 6		8680614811135
NST07	NOSTRIL RETAINER NO: 7		8680614811142
NST08	NOSTRIL RETAINER NO: 8		8680614811159
NST09	NOSTRIL RETAINER NO: 9		8680614811166
NST10	NOSTRIL RETAINER NO: 10		8680614811173
NST11	NOSTRIL RETAINER NO: 11		8680614811180
NST12	NOSTRIL RETAINER NO: 12		8680614811197
NST13	NOSTRIL RETAINER NO: 13		8680614811203